

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 1 of 105

Pre Execution Signatures

Authored by :	Sneha S 00006387 Author	Date:	24 Aug 2026 at 12:54 PM
Reviewed By :	Santhosh Subramanyam 00002803 Reviewer	Date:	24 Aug 2026 at 01:17 PM
Approved By :	Swathi T 00001272 Approver	Date:	24 Aug 2026 at 03:34 PM

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 2 of 105

1. Purpose

To verify the below user requirements:

FR_PROD_COMP_002: The system shall allow a user with the appropriate permission to create product comparisons (except in cases where the product has been nullified) from the Compare Products screen.

FR_PROD_COMP_013: The system shall allow a user to view a graphical IDMP Sections summary in Compare Products, and hover over each segment to see section-level comparison details and counts.

FR_PROD_COMP_014: The system shall allow a user to delete a comparison report from the Compare Products screen, subject to having the required Delete permission.

FR_PROD_COMP_016: The system shall allow a user to filter products in Compare Products using a multi select IDMP Section filter and an associated IDMP Data Quality filter.

Record the following

Product Name	SPORIFY
Product Version	2026.2.0.0
Test Script Approval Status	Approved
Execution Date	01 Sep 2026 at 07:09 PM
Execution Status	Passed
Execution Approval Status	Approved
Run Name[Run ID]	Run_8-27_15-59-44[3504]

2. Pre-Requisites

1. Pre-Validation plan is approved.
2. Tester to complete the required trainings on the application/product/system/ validation testing requirements.
3. ALM Access with respect to project shall be provided.
4. OQ shall be reviewed and approved for execution.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 3 of 105

5. IQ script is executed successfully.
6. MAH configurations are available for the tenant.
7. PMS Synchronisation is completed and target products are available.

3. Test Conditions

Not Applicable

4. Special Instructions

1. The tester shall capture, where possible, objective evidences of the input, the results such as change in status on the status bar, and System generated values, message box, error messages etc., with the System Date and Time stamps of the execution.
2. Steps marked as *PS* shall capture screenshots relevant to that step along with System Date and Time stamps of the execution.
3. The executor shall enter the 'Actual Results' column as "AE" which is an abbreviated for "As Expected" provided the result matches the expected results.
4. 'NA' in the Test Data Column represents that there is no Test Data applicable for that step.

5. Execution Environment


Validation

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 4 of 105


6. Test Steps

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 1	Pre-requisite: 1.User is available with Dashboard and Products domain (full) permission. 2.Product Source List and Source Products are available. 3.Target Products are available *PS*	1.User Name: <pvproduser01> with permissions to Dashboard and Products domain (full) 2.List Name: <SourceList1>	All Pre-requisites are met.	1.User Name: prevaluser3 with permissions to Dashboard and Products domain (full) 2.List Name: Preval-2 AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:03 PM	Passed

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 5 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
 The screenshot displays the SPORIFY web application interface. A modal dialog titled 'Permissions - prevauser3' is open, showing a hierarchical tree of permissions. The tree includes categories such as Pages, Products, Source Data, Source Elements, Source List, Target Products, PMS, Referentials, Substances, and Support. Each category has a list of sub-permissions with checkboxes. The 'SAVE' button at the bottom right of the dialog is highlighted. The background shows the 'Users' management page with a table of users and a sidebar with navigation options.					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 8 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 9 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 2	Launch the SPORIFY application and enter the details as mentioned in test data column. Click on LOGIN button. *PS*	URL: <https://preval.sporify.eu/Account/Login> Tenancy Name: <preval> User name or email: <pvproduser01> Password: <Enter valid password of the user>	The SPORIFY application homepage is displayed.	URL: https://preval.sporify.eu/Account/Login Tenancy Name: AGCS User name or email: prevaluser3 AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:05 PM	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 10 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 11 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
	 The screenshot displays the SPORIFY dashboard interface. At the top, there's a navigation bar with the SPORIFY logo and a user profile. The main content area is divided into several sections: MAPPING STATUS OVERVIEW: A bar chart showing mapping progress across different categories. Substances: A red box indicating '6-SRS Available'. Products: A green box indicating 'SPOR Products 355'. Organisations: A blue box indicating 'Organisations Available'. Referentials: A dark blue box indicating 'Referentials Available'. RESOLUTION STATUS OVERVIEW: Four donut charts showing the status of Substances (Not Set 100%), Products (Auto sync status 96%), Organisations (Not Set 100%), and Referentials (Not Set 100%). TARGET SYSTEM UPDATES: A line chart showing updates over time from Sep-2025 to Jul-2026. SPORIFY ACTIVITY: A section at the bottom showing activity for All and Products. TARGET SYSTEM ACTIVITY: A line chart showing activity over time, including Favourites. The bottom of the screenshot shows a Windows taskbar with the date 5:05 PM 8/27/2026.				

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 12 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
 The screenshot displays the SPORIFY dashboard interface. At the top, there's a navigation bar with the SPORIFY logo and a search bar. Below this, a sidebar menu lists 'Dashboard', 'Products', 'User Guides', and 'Policies'. The main content area is divided into several sections: RESOLUTION STATUS OVERVIEW: Four donut charts showing the status of different entities: Substances (Not Set, 100%), Products (Auto sync status, 96%), Organisations (Not Set, 100%), and Referentials (Not Set, 100%). TARGET SYSTEM UPDATES: A line graph showing updates over time from Sep-2025 to Jul-2026. The graph shows a sharp initial peak followed by a decline and then a series of smaller peaks. SPORIFY ACTIVITY: A list of recent activities including 'PMS Lists Uploaded', 'Products Uploaded', and 'Product Status: Multiple Choice' with associated dates and counts. TARGET SYSTEM ACTIVITY: A list of recent system updates, such as 'Updated PMS Padcev 30 mg powder for concentrate for solution for infusion' and 'Updated PMS Xtandi - 40 mg soft capsules', with their respective dates. The bottom of the dashboard shows the current time zone as (UTC+00:00) Dublin, Edinburgh, Lisbon, London, and the current time as 5:05 PM on 8/27/2026.					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 13 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 3	Expand Products and click on Compare Products. *PS*	NA	New tab is opened and Compare Products screen is displayed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:07 PM	Passed



The screenshot displays the SPORIFY web application interface. The top navigation bar includes the SPORIFY logo and a user profile icon. The main dashboard area is divided into several sections: a left sidebar with navigation links (Dashboard, Products, Overview, Source Lists, Target Products, Compare Products, User Guides, Policies), a central 'Dashboard Statistics and reports' section, and a right sidebar with quick links to Substances, Products, Organisations, and Referentials. The central section contains a 'MAPPING STATUS OVERVIEW' bar chart, a 'RESOLUTION STATUS OVERVIEW' with four donut charts, and two activity charts at the bottom. The Windows taskbar at the bottom shows the system clock as 5:06 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 15 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 4	FR_PROD_COMP_002: Click on the Compare Products button. *PS*	NA	Product Comparison Pop up is displayed with Cancel and Compare buttons.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:08 PM	Passed




	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 16 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 5	Verify fields in the Product Comparison Pop up is displayed as mentioned in the Test data. *PS*	1. Select Primary Product: Editable 2. Select Product: Editable 3. Select Secondary Product: Editable 4. Select Product: Editable	Fields in the Product Comparison Pop up are displayed as mentioned in the Test data.	1. Select Primary Product: Editable 2. Select Product: Editable 3. Select Secondary Product: Editable 4. Select Product: Editable AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:10 PM	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 17 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 18 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 19 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 20 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 21 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 6	Select Primary Product and Product values as mentioned in test data column and verify Product Version field is not displayed for Primary Product. *PS*	1. Select Primary Product: <SourceList1> 2. Select Product: <srcproduct1>	Primary Product details are selected and on selection, Select Product Version field is not displayed for Primary Product.	1. Select Primary Product:Preval-2 2. Select Product: Xospata 40 mg film-coated tablets AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:13 PM	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 22 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 23 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 7	Select Secondary Product and Product values as mentioned in test data column and verify Product Version field is not displayed for Secondary Product. *PS*	1. Select Secondary Product: <SourceList1> 2. Select Product: <srcproduct2>	Secondary Product details are selected and on selection Select Product Version field is not displayed for Secondary Product.	1. Select Secondary Product:Preval-2 2. Select Product: Betmiga 50 mg prolonged-release tablets AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:13 PM	Passed

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 24 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 25 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 8	Click on the Compare button. *PS*	NA	Comparison is initiated.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:15 PM	Passed




The screenshot shows the SPORIFY web application interface. The browser address bar displays the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The user is logged in as AGCS\prevaluser3. The left sidebar contains a menu with the following items: Products (selected), Overview, Source Lists, Target Products, Compare Products, User Guides, and Policies. The main content area is titled 'Compare Products' and includes a search bar and a table of product comparisons. The table has the following columns: Actions, Primary Product, Secondary Product, Export, IDMP Sections, Insights, and Created By. The table contains four rows of data, each representing a product comparison. The first row shows a comparison between Nyracta (SRD) 4 mg - Film-coated tablets and Advagraf 0.5 mg prolonged-release tablets. The second row shows a comparison between Nyracta (SRD) 4 mg - Film-coated tablets and Nyracta (SRD) 4 mg - Film-coated tablets. The third row shows a comparison between Losec Control 20 mg gastro-resistant tablets and Losec Control 20 mg gastro-resistant tablets. The fourth row shows a comparison between Losec Control 20 mg gastro-resistant tablets and prd2. A 'Marketing Authorisation' tooltip is visible over the table. At the bottom of the main content area, there is a blue button that says 'Comparison in Progress'.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 26 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status																																										
SPORIFY SPORIFY + ← ↺ 🏠 https://sporify-preval.sporify.eu/app/main/products/compare-products 🔍 ☆ ☆ 👤 ⋮ Chat SPORIFY AGCS\prevaluser3 Products Overview Source Lists Target Products Compare Products User Guides Policies Compare Products Compare a Primary and Secondary Product Compare Products Q Search keyword × 🔍 <table><tr><th>Actions</th><th>Primary Product 🔼</th><th>Secondary Product 🔼</th><th>Export</th><th>IDMP Sections</th><th>Insights</th><th>Created</th></tr><tr><td></td><td>Xospata 40 mg film-coated tablets</td><td>Botmiga 50 mg prolonged-release tablets</td><td></td><td>Processing</td><td></td><td>pre</td></tr><tr><td>⋮</td><td>Nyracta (SRD) 4 mg - Film-coated tat</td><td>Advagraf 0.5 mg prolonged-release hard capsuli</td><td></td><td>Comparison Unsucce</td><td></td><td>sa</td></tr><tr><td>⋮</td><td>Nyracta (SRD) 4 mg - Film-coated tat</td><td>Nyracta (SRD) 4 mg - Film-coated tablet</td><td></td><td>Comparison Unsucce</td><td></td><td>ad</td></tr><tr><td>⋮</td><td>Lossec Control 20 mg gastro-resistant</td><td>Lossec Control 20 mg gastro-resistant tablets</td><td>📄</td><td> </td><td>[.d]</td><td>ad</td></tr><tr><td>⋮</td><td>Lossec Control 20 mg gastro-resistant</td><td>prd2</td><td>📄</td><td> </td><td>[.d]</td><td>ad</td></tr></table> Total: 5 << < 1 > >> 50 SPORIFY 2026.2.0.0 20260825152117						Actions	Primary Product 🔼	Secondary Product 🔼	Export	IDMP Sections	Insights	Created		Xospata 40 mg film-coated tablets	Botmiga 50 mg prolonged-release tablets		Processing		pre	⋮	Nyracta (SRD) 4 mg - Film-coated tat	Advagraf 0.5 mg prolonged-release hard capsuli		Comparison Unsucce		sa	⋮	Nyracta (SRD) 4 mg - Film-coated tat	Nyracta (SRD) 4 mg - Film-coated tablet		Comparison Unsucce		ad	⋮	Lossec Control 20 mg gastro-resistant	Lossec Control 20 mg gastro-resistant tablets	📄		[.d]	ad	⋮	Lossec Control 20 mg gastro-resistant	prd2	📄		[.d]	ad
Actions	Primary Product 🔼	Secondary Product 🔼	Export	IDMP Sections	Insights	Created																																									
	Xospata 40 mg film-coated tablets	Botmiga 50 mg prolonged-release tablets		Processing		pre																																									
⋮	Nyracta (SRD) 4 mg - Film-coated tat	Advagraf 0.5 mg prolonged-release hard capsuli		Comparison Unsucce		sa																																									
⋮	Nyracta (SRD) 4 mg - Film-coated tat	Nyracta (SRD) 4 mg - Film-coated tablet		Comparison Unsucce		ad																																									
⋮	Lossec Control 20 mg gastro-resistant	Lossec Control 20 mg gastro-resistant tablets	📄		[.d]	ad																																									
⋮	Lossec Control 20 mg gastro-resistant	prd2	📄		[.d]	ad																																									
🪟 🔍 Search 📀 🌐 🔍 📧 📅 📁 🔧 👤 🟢 ⤴ 📶 🔊 🔌 5:14 PM 8/27/2026																																															

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 27 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 9	Verify the above initiated Compare Product record is displayed on top of the table with Comparison status as 'Processing' in IDMP Sections column. *PS*	NA	The above initiated Compare Product record is displayed on top of the table with Comparison status as 'Processing' in IDMP Sections column.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:15 PM	Passed



The screenshot shows the SPORIFY web application interface. The browser address bar displays the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The user is logged in as AGCS\prevaluser3. The interface features a sidebar with navigation options: Products, Source Lists, Target Products, Compare Products (selected), User Guides, and Policies. The main content area is titled 'Compare Products' and includes a search bar and a table of product comparisons. The table has columns for Actions, Primary Product, Secondary Product, Export, IDMP Sections, Insights, and Created. The table lists several products, including Xospata, Nyraeta, and Losec Control, with their respective comparison statuses and IDMP sections. A pagination bar at the bottom indicates 'Total: 5' and shows the current page '1' out of '50'.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 28 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 10	Verify the record is auto refreshed once comparison process is completed. *PS*	NA	Record is auto refreshed once comparison process is completed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:16 PM	Passed



The screenshot shows the SPORIFY web application interface. The browser address bar displays the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The application header includes the SPORIFY logo and the user profile 'AGCS\prevaluser3'. The left sidebar contains navigation links: 'Products', 'Source Lists', 'Target Products', 'Compare Products' (highlighted), 'User Guides', and 'Policies'. The main content area is titled 'Compare Products' and features a search bar and a table of product comparisons. The table has columns: 'Actions', 'Primary Product', 'Secondary Product', 'Export', 'IDMP Sections', 'Insights', and 'Created'. The table lists five comparisons, including Xospata, Nyracta, and Losec Control. A pagination bar at the bottom indicates 'Total: 5' and shows the current page '1' out of '50'.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 29 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 11	Verify IDMP Sections Details are displayed after the comparison process is Completed. *PS*	NA	IDMP Sections Details are displayed after comparison process is Completed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:17 PM	Passed



The screenshot displays the SPORIFY web application interface. The browser address bar shows the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The user is logged in as AGCS\prevaluser3.

The application interface includes a sidebar with the following navigation items:

- Products
 - Overview
 - Source Lists
 - Target Products
 - Compare Products (Active)
- User Guides
- Policies

The main content area is titled "Compare Products" and contains a search bar and a table of product comparisons. The table has the following columns: Actions, Primary Product, Secondary Product, Export, IDMP Sections, and Insights.


Actions	Primary Product	Secondary Product	Export	IDMP Sections	Insights
⋮	Xospata 40 mg film-coated tablets	Botmiga 50 mg prolonged-release tablets	Download		
⋮	Nyracta (SRD) 4 mg - Film-coated tat	Advagraf 0.5 mg prolonged-release hard capsule		Comparison Unsuccessful ⚠	
⋮	Nyracta (SRD) 4 mg - Film-coated tat	Nyracta (SRD) 4 mg - Film-coated tablet		Comparison Unsuccessful ⚠	
⋮	Losec Control 20 mg gastro-resistant	Losec Control 20 mg gastro-resistant tablets	Download		
⋮	Losec Control 20 mg gastro-resistant	prd2	Download		

At the bottom of the table, it indicates "Total: 5" with pagination controls showing page 1 of 50.

The Windows taskbar at the bottom shows the system clock as 5:17 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 30 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 12	Enter the data in the Search field as mentioned in test data column and press [Enter] key. *PS*	<srcproduct1>	Compare products record satisfying the search criteria is displayed and highlighted in yellow.	Xospata 40 mg film-coated tablets AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:20 PM	Passed



The screenshot displays the SPORIFY web application interface. The main content area is titled 'Compare Products' and shows a search for 'Xospata 40 mg film-coated tablets'. The results are presented in a table with columns for 'Actions', 'Primary Product', 'Secondary Product', 'Export', 'IDMP Sections', and 'Insights'. The primary product is 'Xospata 40 mg film-coated tablets' and the secondary product is 'Botmiga 50 mg prolonged-release tablets'. The interface also includes a sidebar with navigation options like 'Products', 'Source Lists', and 'Target Products'. The bottom of the screenshot shows a Windows taskbar with the date 8/27/2026 and time 5:20 PM.


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 31 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 13	FR_PROD_COMP_013: Hover on the IDMP Sections for each segment to see section-level comparison details and counts. *PS*	<Section: Medicinal Products & Attachments Missing: 10 Differences: 1 No Inconsistencies: 2>	Different IDMP sections and it respective counts are displayed.	Section: Medicinal Products & Attachments Missing: 18 Differences: 1 No Inconsistencies: 19 Section: Marketing Authorisation Missing: 22 Section: Clinical Particulars No Inconsistencies: 3 Section: Ingredient No Inconsistencies: 26 Section: Pharmaceutical Product No Inconsistencies: 7 Section: Packaged Medicinal Product Missing: 1 No Inconsistencies: 16 AE	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 32 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
				Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:29 PM	


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 33 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 35 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 36 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 37 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status												
SPORIFY AGCS\prevaluser3 Products Overview Source Lists Target Products Compare Products User Guides Policies Compare Products Compare a Primary and Secondary Product Xospata 40 mg film-coated tablets <table><thead><tr><th>Actions</th><th>Primary Product </th><th>Secondary Product </th><th>Export</th><th>IDMP Sections</th><th>Insight</th></tr></thead><tbody><tr><td></td><td>Xospata 40 mg film-coated tablets</td><td>Botmiga 50 mg prolonged-release tablets</td><td></td><td></td><td></td></tr></tbody></table> Marketing Authorisation Total: 1						Actions	Primary Product	Secondary Product	Export	IDMP Sections	Insight		Xospata 40 mg film-coated tablets	Botmiga 50 mg prolonged-release tablets			
Actions	Primary Product	Secondary Product	Export	IDMP Sections	Insight												
	Xospata 40 mg film-coated tablets	Botmiga 50 mg prolonged-release tablets															


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 41 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
	The screenshot displays the SPORIFY application interface. The top navigation bar includes the SPORIFY logo and a user profile 'AGCS\prevaluser3'. The left sidebar contains a 'Products' menu with sub-items: Overview, Source Lists, Target Products, Compare Products (highlighted), User Guides, and Policies. The main content area is titled 'Compare Products' and shows a search for 'Xospata 40 mg film-coated tablets'. Below the search bar is a table with columns: Actions, Primary Product, Secondary Product, Export, IDMP Sections, and Insights. The table lists two products: 'Xospata 40 mg film-coated tablets' and 'Botmiga 50 mg prolonged-release tablets'. The IDMP Sections column shows a bar chart with various colored segments. A tooltip labeled 'Ingredient' is visible over the bar chart. The bottom of the screen shows the Windows taskbar with the date and time as 5:27 PM on 8/27/2026.				


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 42 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 43 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 44 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 45 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 46 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 48 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 14	Click on the Insights icon and verify. *PS*	NA	Data Quality Insights Pop up is displayed and all IDMP Sections details with comparison counts are displayed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:30 PM	Passed




The screenshot displays the SPORIFY web application interface. The browser address bar shows the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products?search=Xospata%2040%20mg%20film-coated%20tablets&searchTyp...>. The application header shows the user is logged in as AGCS\prevaluser3. The left sidebar contains navigation links: Products, Overview, Source Lists, Target Products, Compare Products (active), User Guides, and Policies. The main content area is titled 'Compare Products' and includes a search bar with the text 'Xospata 40 mg film-coated tablets'. Below the search bar, a table displays the comparison results. The table has columns for 'Secondary Product 1', 'Export', 'IDMP Sections', 'Insights', and 'Created By'. The first row shows 'Betmiga 50 mg prolonged-release tablets' with a bar chart in the 'IDMP Sections' column and a 'Data Quality Insights' pop-up. The pop-up shows a bar chart and a table with the following data:


Section	Count
Primary	1
Secondary	1
Comparison	1
Total	3

The bottom of the screen shows a Windows taskbar with the time 5:29 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 49 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 50 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 51 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 15	Close the Data Quality Insights Popup.	NA	Data Quality Insights popup is closed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:30 PM	Passed
Step 16	Click on the Export Comparison Report button and verify. *PS*	NA	Comparison Report is Exported to Excel.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:31 PM	Passed

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 52 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
	SPORIFY AGCS\prevaluser3 Products Overview Source Lists Target Products Compare Products User Guides Policies Compare Products Compare a Primary and Secondary Product Xospata 40 mg film-coated tablets Secondary Product T1 Export IDMP Sections Insights Created By ted tablets Betmiga 50 mg prolonged-release tablets Export Comparison Report Total: 1 << < 1 > >> 50				

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 53 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 54 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 17	Open the exported comparison report and verify. *PS*	NA	Exported comparison report contains details of Primary and Secondary products and the summary of comparison details.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:35 PM	Passed

PROTECTED VIEW Be careful—files from the Internet can contain viruses. Unless you need to edit, it's safer to stay in Protected View. Enable Editing					
B2					38
	A	B	C	D	E
1	Resource	Number of elements	Number of different values between {PMS} vs {RIM}	Number of missing values in {PMS}	Number of missing values in {RIM}
2	Medicinal Product	38	1	9	9
3	Marketing authorisation	22	0	11	11
4	Clinical particulars	3	0	0	0
5	Packaged medicinal product	17	0	1	0
6	Ingredient	26	0	0	0
7	Pharmaceutical product	7	0	0	0
8	Total	113	1	21	20
9					
10					
11	PMS filename	68_1653_1_4.4.json			
12	PMS data refresh	27-08-2026			
13	RIM filename	68_1655_1_4.4.json			
14	RIM data refresh	27-08-2026			
15	Report creation	27-08-2026			
16	Timezone	UTC			
17	Username	sporifyadmin			
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					
28					
Comparison Summary Comparison Report					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 55 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
 The screenshot displays a 'PROTECTED VIEW' window with a warning: 'Be careful—files from the Internet can contain viruses. Unless you need to edit, it's safer to stay in Protected View.' Below the warning is a table comparing two medicinal products. The table has columns for 'Resource', 'IG Reference', 'Element Name', 'PMS', and 'RIM'. The 'PMS' column contains values for 'Medicinal Product Definition/6000000001885' and 'Medicinal Product Definition/6000000000155'. The 'RIM' column contains values for 'Medicinal Product Definition/6000000000155'. The table lists various attributes such as 'Product Management Service Identifier (PMS ID)', 'Domain', '(Authorised) Pharmaceutical Form', 'Combined Pharmaceutical Dose Form', 'Legal Status of Supply', 'Additional Monitoring Indicator', 'PaediatricUseIndicator', 'Product classification', 'xEMPD product type information', 'Legal basis', 'ATC code', 'Genetically Modified Organisms (GMOs)', 'Medicinal product name', 'Full name', 'Country/Language', 'Country', 'Language', '(Medicinal product name) name part(s)', 'Name part type', 'Name part text', and '(Medicinal product name) name part(s)'. The 'Actual Results' column shows 'Missing' for several attributes, including 'Medicinal product name', 'Full name', 'Country/Language', 'Country', 'Language', '(Medicinal product name) name part(s)', 'Name part type', 'Name part text', and '(Medicinal product name) name part(s)'. The 'Expected Results' column shows the corresponding values for the two products.					

[Click Here](#)

Step 18	Click on Clear filter icon.	NA	Search criteria is cleared and all compared products are listed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:36 PM	Passed
---------	-----------------------------	----	--	--	--------

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 56 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 19	Click on the Compare Products button. *PS*	NA	Product Comparison Pop up is displayed with Cancel and Compare buttons.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:36 PM	Passed



	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 57 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 20	Select Primary Product and Product values as mentioned in test data column and verify Product Version field is not displayed for Primary Product. *PS*	1. Select Primary Product: <SourceList1> 2. Select Product: <srcprod1>	Primary Product details are selected and on selection Select Product Version field is not displayed for Primary Product.	1. Select Primary Product:Preval-2 2. Select Product:Betmiga 25 mg prolonged-release tablets AE Chandrashekar Naguvannahalli Shankaregowda 27 Aug 2026 at 05:40 PM	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 58 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 59 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 21	Select Secondary Product and Product values as mentioned in test data column and verify Product Version field is displayed for Secondary Product. *PS*	1. Select Secondary Product: PMS 2. Select Product: <pmsprod1>	Secondary Product details are selected and on selection Select Product Version field is displayed for Secondary Product.	1. Select Secondary Product: PMS 2. Select Product: Xtandi - 40 mg film-coated tablets AE Chandrashekar Naguvannahalli Shankaregowda 27 Aug 2026 at 05:42 PM	Passed

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 60 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 61 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 22	Select a version for the Secondary Product from the dropdown list. *PS*	Select Product Version: <3>	Product version is selected for secondary product.	Select Product Version: 248 AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:43 PM	Passed



	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 62 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 23	Click on Swap icon available between Primary and Secondary product section. *PS*	NA	Primary and Secondary Products are swapped each other along with product and version details.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:43 PM	Passed




The screenshot displays the SPORIFY web application interface. A 'Product Comparison' dialog box is open, allowing users to compare two products. The dialog box contains the following fields:

- Select Primary Product:** PMS
- Select Secondary Product:** Preval-2
- Select Product:** Xtandi - 40 mg film-coated tablets
- Select Product:** Betmiga 25 mg prolonged-release tablets
- Select Product Version:** 248

At the bottom of the dialog box, there are 'Cancel' and 'Compare' buttons. The background shows the SPORIFY application with a sidebar menu and a table of products. The table has columns for 'Insights' and 'Created By'. The 'Created By' column lists users: prevaluser3, santhosh, adminSPORIFY, and adminSPORIFY. The bottom of the screen shows the Windows taskbar with the time 5:43 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 63 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 24	Click on the Compare button. *PS*	NA	Comparison is initiated.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:44 PM	Passed




The screenshot shows the SPORIFY web application interface. The browser address bar displays the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The user is logged in as AGCS\prevaluser3. The main content area is titled 'Compare Products' and contains a table with the following data:

Secondary Product	Export	IDMP Sections	Insights	Created By
ted tablets	Export	Processing		prevaluser3
ted tablets	Export	Comparison Unsuccessful		prevaluser3
n-coated tat	Export	Comparison Unsuccessful		santhosh
n-coated tat	Export	Comparison Unsuccessful		adminSPORIFY
stro-resistant	Export	Comparison Unsuccessful		adminSPORIFY
stro-resistant	Export	Comparison Unsuccessful		adminSPORIFY

The interface also includes a sidebar with navigation options: Products, Source Lists, Target Products, Compare Products, User Guides, and Policies. A 'Comparison In Progress' button is visible at the bottom right.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 64 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 25	Verify the above initiated Compare Product record is displayed on top of the table with Comparison status as 'Processing' in IDMP Sections column. *PS*	NA	The above initiated Compare Product record is displayed on top of the table with Comparison status as 'Processing' in IDMP Sections column.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:45 PM	Passed



The screenshot displays the SPORIFY web application interface. The top navigation bar shows the SPORIFY logo and the user 'AGCS\prevaluser3'. The left sidebar contains a menu with options: Products, Source Lists, Target Products, Compare Products (highlighted), User Guides, and Policies. The main content area is titled 'Compare Products' and features a search bar and a table of product comparisons. The table has columns for Actions, Primary Product, Secondary Product, Export, IDMP Sections, and Insights. The table lists several product comparisons, including Xtandi, Xospata, Nyracta, and Losec Control. The 'IDMP Sections' column shows a progress bar and status indicators. The bottom of the screen shows a Windows taskbar with the date and time as 5:44 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 65 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 26	Verify the record is auto refreshed once comparison process is completed. *PS*	NA	Record is auto refreshed once comparison process is completed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:46 PM	Passed



The screenshot displays the SPORIFY web application interface. The browser window shows the URL <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The application has a sidebar with navigation options: Products, Source Lists, Target Products, Compare Products (highlighted), User Guides, and Policies. The main content area is titled 'Compare Products' and includes a search bar and a table of product comparisons. The table has columns for Actions, Primary Product, Secondary Product, Export, IDMP Sections, and Insights. The table lists several product comparisons, including Xtandi, Xospata, Nyracta, and Losec Control. A green notification bar at the bottom right indicates 'Comparison Complete'. The browser's taskbar at the bottom shows the time as 5:45 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 66 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 27	Verify IDMP Sections Details are displayed after the comparison process is Completed. *PS*	NA	IDMP Sections Details are displayed after the comparison process is Completed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:46 PM	Passed



The screenshot displays the SPORIFY web application interface. The browser address bar shows the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The user is logged in as AGCS\prevaluser3. The interface features a sidebar with navigation options: Products, Source Lists, Target Products, Compare Products (highlighted), User Guides, and Policies. The main content area is titled 'Compare Products' and includes a search bar and a table of product comparisons. The table has columns for Actions, Primary Product, Secondary Product, Export, IDMP Sections, and Insights. The table lists several products, including Xtandi, Xospata, Nyrocta, and Losec Control, showing their respective IDMP section details and comparison status. The bottom of the screenshot shows the Windows taskbar with the system clock at 5:46 PM on 8/27/2026.


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 67 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 28	FR_PROD_COMP_013: Hover on the IDMP Sections for each segment to see section-level comparison details and counts. *PS*	<Section: Medicinal Products & Attachments Missing: 10 Differences: 1 No Inconsistencies: 2>	Different IDMP sections and its respective counts are displayed.	Section: Medicinal Products & Attachments Missing: 46 Differences: 4 No Inconsistencies: 9 Section: Manufacturer Missing: 96 Section: Marketing Authorisation Missing: 25 Section: Clinical Particulars Missing: 4 Section: Ingredient Missing: 82 Section: Pharmaceutical Product Missing: 20 Section: Packaged Medicinal Product Missing: 59 AE	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 68 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
				Chandrashekar Naguvannahalli Shankaregowda 27 Aug 2026 at 05:55 PM	


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 70 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 71 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 72 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 79 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 81 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 85 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 29	Click on the Insights icon and verify. *PS*	NA	Data Quality Insights Pop up is displayed and all IDMP Sections details with comparison counts are displayed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:56 PM	Passed




The screenshot displays the SPORIFY web application interface. The browser address bar shows the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The application header includes the SPORIFY logo and the user profile 'AGCS\prevaluser3'. The left sidebar contains a menu with the following items: Products, Overview, Source Lists, Target Products, Compare Products (highlighted), User Guides, and Policies. The main content area is titled 'Compare Products' and features a search bar and a table of products. The table has columns for 'Secondary Product', 'Export', 'IDMP Sections', 'Insights', and 'Created By'. The products listed include 'n-coated tablets', 'm-coated tablets', 'g - Film-coated tab', and 'ng gastro-resistant'. A 'Data Quality Insights' pop-up is displayed over the table. The bottom of the screenshot shows a Windows taskbar with the date and time as 5:55 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 86 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 87 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 88 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 30	Close the Data Quality Insights Popup.	NA	Data Quality Insights popup is closed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:56 PM	Passed
Step 31	Click on the Export Comparison Report button and verify. *PS*	NA	Comparison Report is Exported to Excel.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:56 PM	Passed

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 89 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 90 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 91 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 32	Open the exported comparison report and verify. *PS*	NA	Exported comparison report contains details of Primary and Secondary products and the summary of comparison details.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 05:58 PM	Passed


PROTECTED VIEW Be careful—files from the Internet can contain viruses. Unless you need to edit, it's safer to stay in Protected View. Enable Editing

A	B	C	D	E
Resource	Number of elements	Number of different values between {PMS} vs {RIM}	Number of missing values in {PMS}	Number of missing values in {RIM}
Medicinal Product	59	4	30	16
Manufacturer	96	0	96	0
Marketing authorisation	25	0	14	11
Clinical particulars	4	0	1	3
Packaged medicinal product	59	0	43	16
Ingredient	82	0	56	26
Pharmaceutical product	20	0	13	7
Total	345	4	253	79
PMS filename	68_1652_1_4.4.json			
PMS data refresh	27-08-2026			
RIM filename	600000000083_248_4.4.json			
RIM data refresh	27-08-2026			
Report creation	27-08-2026			
Timezone	UTC			
Username	sporifyadmin			

Comparison Summary Comparison Report

5:57 PM 8/27/2026

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 92 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

[Click Here](#)

Step 33	FR_PROD_COMP_014: Click on the Actions icon (3 dots) beside the recently compared record and click on Delete Comparison option. *PS*	Primary Product: <srcprod1> Secondary Product: <pmsprod1>	A message "Are You Sure? The comparison between <primary_product> and <secondary_product> will be deleted." is displayed.	Primary Product: Xtandi - 40 mg film-coated tablets Secondary Product: Betmiga 25 mg prolonged-release tablets AE	Passed
---------	--	--	---	---	--------

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 93 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
				Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 06:00 PM	

SPORIFY

Products

Overview

Source Lists

Target Products

Compare Products

User Guides

Policies

AGCS\prevaluser3

Compare Products

Compare a Primary and Secondary Product

Compare Products

Search keyword

Actions	Primary Product	Secondary Product	Export	IDMP Sections	Insights
	Xtandi - 40 mg film-coated tablets	Botmiga 25 mg prolonged-release tablets			
Delete Comparison		Botmiga 50 mg prolonged-release tablets			
	Nyracta (SRD) 4 mg - Film-coated tablet	Advagraf 0.5 mg prolonged-release hard capsules		Comparison Unsuccessful	
	Nyracta (SRD) 4 mg - Film-coated tablet	Nyracta (SRD) 4 mg - Film-coated tablet		Comparison Unsuccessful	
	Losec Control 20 mg gastro-resistant	Losec Control 20 mg gastro-resistant tablets			
	Losec Control 20 mg gastro-resistant	prd2			

Total: 6


javascript;

SPORIFY 2026.2.0.0 20260825152117

Search


5:59 PM 8/27/2026

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 94 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 95 of 105


Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 34	Click on Yes button. *PS*	NA	Selected Comparison record is deleted from the list.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 06:00 PM	Passed



The screenshot displays the SPORIFY web application interface. The browser window shows the URL `https://sporify-preval.sporify.eu/app/main/products/compare-products`. The application header includes the SPORIFY logo and the user profile `AGCS\prevaluser3`. The left sidebar contains navigation links: **Products** (Overview, Source Lists, Target Products, Compare Products), **User Guides**, and **Policies**. The main content area is titled **Compare Products** and features a search bar and a table of product comparisons. The table has columns for **Actions**, **Primary Product**, **Secondary Product**, **Export**, **IDMP Sections**, and **Insight**. The table lists five comparisons, including products like Xospata, Nyracta (SRD), and Losec Control. A green notification banner at the bottom right indicates **Comparison Deletion Successful**. The footer of the application shows the version `SPORIFY 2026.2.0.0 20260825152117`. The Windows taskbar at the bottom shows the time as 6:00 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 96 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 35	FR_PROD_COMP_016: Click on the Advanced Search icon. *PS*	NA	Advanced Search flyout menu is displayed with only IDMP Sections filter.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 06:01 PM	Passed



The screenshot displays the SPORIFY web application interface. The browser address bar shows the URL: <https://sporify-preval.sporify.eu/app/main/products/compare-products>. The application has a sidebar on the left with the following menu items: Products (selected), Overview, Source Lists, Target Products, Compare Products, User Guides, and Policies. The main content area is titled 'Compare Products' and includes a search bar labeled 'Search keyword'. Below the search bar, there is a table with columns 'Actions', 'Primary Product', and 'Secondary Product'. The table contains the following data rows:


Actions	Primary Product	Secondary Product
	Xospata 40 mg film-coated tablets	Botmiga 50 mg prolonged-release tablets
	Nyracta (SRD) 4 mg - Film-coated tat	Advagraf 0.5 mg prolonged-release hard caps
	Nyracta (SRD) 4 mg - Film-coated tat	Nyracta (SRD) 4 mg - Film-coated tablet
	Losec Control 20 mg gastro-resistant	Losec Control 20 mg gastro-resistant tablets
	Losec Control 20 mg gastro-resistant	prd2

On the right side of the main content area, there are buttons for 'Save View as Global Default' and 'Reset to Default'. Below these buttons, there are tabs for 'Preferences' and 'Advanced Search'. Under 'Advanced Search', there is a section titled 'Domain Filters - Compare Products' with a dropdown menu labeled 'IDMP Sections' and a text input field 'Select IDMP Section(s)'. The bottom of the screenshot shows a Windows taskbar with a search bar, several application icons, and a system clock showing 6:01 PM on 8/27/2026.

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 97 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 36	Enter the data in the IDMP Sections filter as mentioned in test data column. *PS*	IDMP Sections: <Marketing Authorisation, Pharmaceutical Product>	Values are selected for IDMP Sections filter and another filter IDMP Data Quality is displayed.	IDMP Sections: Marketing Authorisation, Pharmaceutical Product AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 06:03 PM	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 98 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 99 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 37	Enter the data in the IDMP Data Quality filter as mentioned in test data column and click on background window. *PS*	IDMP Data Quality: <Missing, No Inconsistencies>	Values are selected for IDMP Data Quality filter and comparison records satisfying the search criteria are displayed.	IDMP Data Quality: Missing, No Inconsistencies AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 06:04 PM	Passed


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 100 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					


	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 102 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
Step 38	Navigate to previous tab, Click on User Account details and click on Logout option. *PS*	NA	The SPORIFY application Log in page is displayed.	AE Chandrashekar Naguvanahalli Shankaregowda 27 Aug 2026 at 06:05 PM	Passed

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 103 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
	 The screenshot displays the SPORIFY dashboard interface. At the top, there's a navigation bar with the SPORIFY logo and a sidebar menu on the left containing options like Dashboard, Products, Overview, Source Lists, Target Products, Compare Products, User Guides, and Policies. The main content area is titled 'Dashboard Statistics and reports' and includes several widgets: a 'MAPPING STATUS OVERVIEW' bar chart, a 'RESOLUTION STATUS OVERVIEW' with four donut charts for Substances (Not Set 100%), Products (Auto sync status 96%), Organisations (Not Set 100%), and Referentials (Not Set 100%), a 'TARGET SYSTEM UPDATES' line chart, and a 'SPORIFY ACTIVITY' section at the bottom. A user profile dropdown menu is visible in the top right corner, showing options like 'Change password', 'Login attempts', 'Change profile picture', 'My settings', and 'Logout'. The bottom of the image shows a Windows taskbar with the time 6:04 PM on 8/27/2026.				

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 104 of 105

Step No	Procedure Steps	Test Data	Expected Results	Actual Results	Test Status
					

	Project	SPORIFY	Document Type	OQ Script
			Document Version	1.0
Document Name	AG_SPOR_2026.2.0.0_OQ_Products_Comparison		Page	Page 105 of 105

7. Post Execution Signatures

Executed By :	Chandrashekar Naguwanahalli Shankaregowda 00005325 Executor	Date:	01 Sep 2026 at 07:09 PM
Reviewed By :	Santhosh Subramanyam 00002803 Reviewer	Date:	02 Sep 2026 at 04:28 PM
Approved By :	Swathi T 00001272 Approver	Date:	03 Sep 2026 at 11:01 AM